S*PURE Application Note - Environment

PAH's in Water On SEClute SPE Column



Procedure using SEClute™ C18-Max SPE tube, 500mg/6mL (Part #5138767)

Sample Treatment: 500mL of water and add 75mL of 2-propanol, spiked with 20uL of PAH mixture at 2000ug/mL of each PAH, which corresponds to 40ug/500mL water, or 80ppb of each PAH.

Conditioning: 2 x 5mL of methanol, then 5mL of 85% water/15% 2-propanol.

Sample Application: Attach 75mL reservoir (Part #210575) with a syringe adapter (Part #210705) and apply entire sample. Dry columns for 30 minutes at full vacuum.

Elution: Elute PAHs with 3 x 1mL CH2Cl2 and adjusted the volume to 1mL by evaporating under a gentle nitrogen stream.

GC Column: AT 5 ms columns, (Part #15807)

Oven: 80°C (1min hold) to 160°C at 25°C/min, to 300°C at 3°C/min, to 325°C (3 min hold) at 25°C/min


Injector: 270°C

FID: 340°C

Split Ratio: 50:1, 14.0psig, helium Linear Velocity: 28cm/sec (at 80°C)

1.0uL injections were done via

7683B autosampler.



- 1. Naphthalene
- 2. Acenaphthylene
- 3. Acenaphthene 4.Fluorene
- 5.Phenanthrene
- 6.Anthracene
- 7.Fluoranthene 8.Pyrene
- 9. Benzo(a)anthracene
- 10.Chrysene
- 11.Benzo(b)fluoranthene 12.Benzo(k)fluoranthene
- 13.Benzo(a)pyrene
- 14.Indeno(1,2,3-cd)pyrene
- 15.Dibenzo(a,h)anthracene
- 16.Benzo(ghi)perylene

www.sstarpure.com